

test-it[®] 4416

One-piece 3 cells 4x4x16 cm³ moulds

New :

**Mould 4x4x16
for withdrawal
measurement**



test-it[®] 4416R

Safety / Quality / Comfort / Efficiency / Economy

- Specimens do not stick to the sides of the mould when ejected
 - Perfectly flat and smooth sides
 - Easy to clean
 - Rustproof and airtight
 - Light (1 kg)
- Usable on jolting tables or vibration tables
 - Ideal for the very fluid products
- Perfect also with resin cements and mortars



test-it[®] 4416

	test-it® 4416 mould	Detachable mould : steel ¹ or plastic ²	Polystyrene mould
Handling	Easy & Safe	Restrictive and bothersome, even dangerous ¹	Fragile and can be bent out of shape
Finish	Perfectly smooth Specimens that do not stick	Product stick to the mould, breaks easily, leaks	Rough, geometric problems
Cleaning	Effortless Cast in one piece, easy to clean with a sponge or compressed air	Difficult and time-consuming : disassembly, assembly, greasing	Disposed of waste
Life	Several years *	Limited : corrosion ¹ , assembly, brittle ²	Disposable : single use only



test-it® spray **

Special emulsion available in spray.
Ensures easy maintenance and ejection from moulds.



* Subject to use in accordance with the instructions provided with **test-it® spray**

** The use of other cleaning or release agents is procribed and may damage your moulds irremediably



MLB – 7, chemin de Saint Jean – 06650 Le Rouret
Tél : +33 951 246 370 – Fax : +33 956 246 370

contact@mlbsas.eu – www.mlbsas.eu